

Data Logger Sets for Operation and Process Qualification

- **EBI 11 Mini Data Logger Sets**

SL 1221: Validation set for washer-disinfectors for endoscopes (see page 50)

SL 3301: Validation set for DAC Universal, small sterilizers and washer-disinfectors (see page 52)

- **EBI 12 Precision Data Logger Sets**

SL 2001: for washer-disinfector (see page 51)

SL 3001: for benchtop sterilizers (see page 51)

SL 3101: for large steam sterilizers (see page 52)

The sets can individually be expanded or compiled by yourself from one or more data loggers (EBI 11, EBI 16 or EBI 12), the appropriate interface and the corresponding software. Talk to us!

SL 1221 Complete Validation Set for Washer-Disinfectors for Endoscopes



For the validation of thermal processes in washer-disinfectors for endoscopes according to ISO 15883.

This reliable system consists of user-friendly mini data loggers and an Endoscope-Dummy that can be placed directly in the processes.



The set contains:

- 3 x EBI 11-T235 Mini temperature data logger, needle length: 25 mm
- 1 x EBI 11-P111 Mini pressure data logger
- Sealing kit EBI 11-Valiset
- EBI IF 300, 4-port Interface
- Winlog.validation, evaluation software
- spypach Endoscope-Dummy "spo-pro" Professional
- Aluminum carrying case

Adequate pH and conductivity meters, see page 55.

Type	Description	Part No.
SL 1221	Endoscope set validation PROFESSIONAL	1250-1221
Spo-pro BASIC	Dummy BASIC without logger	1340-6086
Spo-pro-CLASSIC	Dummy CLASSIC without logger	1340-6087
Spo-pro PROFESSIONAL	Dummy PROFESSIONAL without logger	1340-6088