

SL 2001 Complete Validation Set for washer-disinfectors



For the validation of washer-disinfectors according to ISO 15883. The set can individually be expanded or compiled yourself from one or more data loggers (EBI 12 or EBI 11), the appropriate interface and corresponding software. Talk to us!

The set contains:

- 1 x EBI 12-T220 Temperature data logger
- 2 x EBI 12-T441 Temperature data loggers
- 1 x EBI 12-TP231 Temperature/pressure data logger
- EBI IF 200, 4-port Interface with USB connection and antenna
- TÜV certified software Winlog.validation
- Carrying case „SYSTAINER“



Type	Description	Part No.
SL 2001	Validation set for washer-disinfectors	1250-2001

SL 3001 Complete Validation Set for benchtop steam sterilizers



For the validation of small steam sterilizers according to ISO 17665.

The set can individually be expanded or compiled yourself from one or more data loggers (EBI 11, EBI 16 or EBI 12), the appropriate interface and corresponding software. Talk to us!

The set contains:

- 1 x EBI 12-TP453 Temperature/pressure data logger with AL 101 silicone protection box
- EBI IF 200 4-port Interface with USB connection and antenna
- TÜV certified software Winlog.validation
- Carrying case „SYSTAINER“



Type	Description	Part No.
SL 3001	Validation set for small steam sterilizers	1250-3001